


## Lambda DNA Restriction Map



Taken from A Sourcebook of Biotechnology Activities, published by the National Association of Biology Teachers, 1990. ISBN 0-941212-09-2

NEBuffer	1	2	3	4
% Activity	75	100	50	75

**Source:** An *E. coli* strain that carries the cloned *BamH I* gene from *Bacillus amyloquelaciens* H (ATCC 49763)

**Reaction Buffer:** NEBuffer *BamH I* + BSA  
150 mM NaCl, 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 1 mM dithiothreitol (pH 7.9 @ 25°C). Supplement with 100 µg/ml BSA. Incubate at 37°C.

**Ligation and Recutting:** After 50-fold overdigestion with *BamH I*, > 95% of the DNA fragments can be ligated and recut.

**Concentration:** 20,000 and 100,000 units/ml. Assayed on λ DNA.

**Storage Conditions:** 50 mM KCl, 10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA, and 50% glycerol. Store at -20°C.



### BamH I

#R0136S 10,000 units ..... \$50  
#R0136L 50,000 units ..... \$200

for high (5X) concentration, order #R0136T (10,000 units) or #R0136M (50,000 units)



2000-01 Catalog & Technical Reference

TELEPHONE ORDERING  
1-800-632-5227  
(1-800-NEB-LABS)  
8:00 to 6:00 pm EST

FAX ORDERING  
1-800-632-7440

E-MAIL ORDERING  
orders@neb.com

TECHNICAL ASSISTANCE  
1-800-632-7799  
info@neb.com

INTERNET/ON-LINE ORDERS  
www.neb.com

USA  
New England Biolabs, Inc.  
32 Tozer Road  
Beverly, MA 01915-5599  
USA  
Telephone (978) 927-5054  
Toll Free (USA Orders) 1-800-632-5227  
Toll Free (USA Tech) 1-800-632-7799  
Fax (978) 921-1350  
e-mail: info@neb.com  
www.neb.com

NEBuffer	1	2	3	4
% Activity	100	100	100	100

**Source:** An *E. coli* strain that carries the cloned *EcoR I* gene from *E. coli* RY13 (R.N. Yoshimori)

**Reaction Buffer:** NEBuffer *EcoR I*  
50 mM NaCl, 100 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 0.025% Triton X-100 (pH 7.5 @ 25°C). Incubate at 37°C.

**Ligation and Recutting:** After 100-fold overdigestion with *EcoR I*, > 95% of the DNA fragments can be ligated and recut.

**Concentration:** 20,000 and 100,000 units/ml. Assayed on λ DNA.



### EcoR I

#R0101S 10,000 units ..... \$50  
#R0101L 50,000 units ..... \$200

for high (5X) concentration, order #R0101T (10,000 units) or #R0101M (50,000 units)



**Storage Conditions:** 300 mM NaCl, 10 mM KPO<sub>4</sub> (pH 7.5), 0.1 mM EDTA, 10 mM 2-mercaptoethanol, 0.15% Triton X-100, 200 µg/ml BSA, and 50% glycerol. Store at -20°C.

**Diluent Compatibility:** Diluent C, see page 205.

**Heat Inactivation:** 65°C for 20 minutes.

**Note:** Conditions of low ionic strength, high enzyme concentration, glycerol concentration >5%, or pH >8.0 may result in star activity (see page 213). For performing double digests with *EcoR I*, see page 204.

### Hind III

#R0104S 10,000 units \$50  
#R0104L 50,000 units \$200

for high (5X) concentration, order #R0104T (10,000 units) or #R0104M (50,000 units)



**Source:** An *E. coli* strain that carries the cloned *Hind III* gene from *Haemophilus influenzae* Rd (ATCC 51907)

**Reaction Buffer:** NEBuffer 2  
50 mM NaCl, 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 1 mM dithiothreitol (pH 7.9 @ 25°C). Incubate at 37°C.

**Ligation and Recutting:** After 200-fold overdigestion with *Hind III*, > 95% of the DNA fragments can be ligated and recut.


**Concentration:** 20,000 and 100,000 units/ml. Assayed on λ DNA.

**Storage Conditions:** 250 mM NaCl, 10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM dithiothreitol, 500 µg/ml BSA, and 50% glycerol. Store at -20°C.

NEBuffer	1	2	3	4
% Activity	50	100	10	50

**Diluent Compatibility:** Diluent B, see page 205.

**Heat Inactivation:** 65°C for 20 minutes.



*Hind III* crystals (Ira Schildkraut and Lydia Dorner, New England Biolabs)

