

PCR Master Mix- Thermopol Buffer (NEB)- All primer sets **except Mu**

A Master Mix is made and used to save time, reagents, and maximize reproducibility and reliability of the protocol and hence, the results.

The PCR Master Mix consists of PCR grade **H₂O**, **PCR buffer**, **MgCl₂** (may be included in the PCR buffer or added separately), **dNTPs**, **primers** (The primers may or may not be included in the Master Mix, depending on the total number of reactions and primer sets used.) and **thermostable polymerase** (i.e. Taq).

A Master Mix is made by multiplying each component by the total number of PCR reactions per primer set (including + and negative controls, and duplicates). In addition, an extra reaction should be added for every 10 reactions to correct for volume errors due to repetitive pipetting.

***Example – 25 rxns (25 +2 = 27)**

Vol. (μl)/X rxns.	Vol. (μl)/1 rxn.	Reagent	Stock Conc.	Final Conc.
135	5	Buffer w/MgCl ₂	10X	1X
27	1	dNTPs	10 mM	0.2 mM
27	1	Forward primer	10 μM	0.2 μM
27	1	Reverse primer	10μM	0.2 μM
6.75	0.25	Polymerase	5 units/μl	1.25 units
992.25	36.75	Water		
1215	45	Total volume		
	5	DNA		5 μl/rxn.

45 μl/rxn.

50 μl total

*** Note:** Total PCR reaction volume and the volume of DNA added to a reaction may vary. These volumes are dependent on the protocol.